

August 26, 2020

Len-Tex has performed multiple rounds of testing on our products to assess their durability after repeated exposure to various institutional disinfectants with different chemistries. Please note the methods and disinfectants were tested in a controlled setting, not necessarily meant to simulate any specific scenario or environment. These tests did NOT result in a recommendation for any specific disinfectants, but rather a report of our findings.

In regards to methodology when applying disinfectants, following the manufacturer's instructions is recommended. Take note of the product's application method, contact time, and whether or not a rinse is required. The disinfectants that were tested by Len-Tex that had no effect on the product are listed below:

- **Nixall**
- **Neutral Q**
- **Isopropyl Alcohol** (60% and 70% concentration)
- **Virex**
- **Avistat D**
- **Clorox Wipes** (non-bleach)
- **Clorox Bleach** (50% concentration) + water rinse